NOVA Cer[®] MC

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Safety information: *Please read these instructions thoroughly before using the product. Improper use of the product and failure to observe the information provided may affect quality and reduce the lifetime of the dental prosthesis. Only trained staff (dentists or dental technicians) should be allowed to handle this product.*

– WARNUNG –

Applies to porcelains in powder form: Please open the packaging and use the product only with protective equipment. Wear safety goggles and a dust mask. Dust particles irritate the eyes, mucous membranes and skin, and injure the lungs.

Intended use: NOVACer® for metalceramics (MC) covers modular ceramic materials for veneering the metal framework .

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Alloy:

Precious metal or non-precious metals with a CTE of 13,8 - 14,9 at 25°C - 500°C can be veneered.

CTE > 14,5: Prolonged cooling

CTE < 14,1: Object must be removed rapidly from the firing chamber

Oxidize according to the alloy manufacturer's instructions.

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Applying NOVACer® NE Bonder:

Stir before using! Apply the bonder in a thin layer but mask the framework completely.

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1. Opaque Firing:

Apply the NOVACer[®] MC paste- or powder opaquer on to the metal frame with approx. 75% coverage. The opaque firing should shine.

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2. Opaque Firing:

The second opaque firing is applied to give 100% coverage of frame and should slightly shine.

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Technical advice to NE Bonder:

NOVACer® NE Bonder is used as a layer between the alloy and the ceramic. It blocks escaping metal oxides and neutralize differences in the thermal expansion coefficient.

If Bonder has dried out, do not re-moisten. Close the lid after using.

Technical advice to paste opaque:

- Do not moisten the brush with water
- Avoid any dilution with water
- The consistency can be modified slightly with glaze liquid
- Avoid drying too quickly
- The pre-drying time must be adhered to

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*Build-up of opaque dentin in the neck area with
NOVACer[®] MC opaque dentine or dentine modifier.*

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The labial or occlusal shape is build-up with NOVACer[®] MC dentine in the required shade.

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The dentin will be reduced to create enough space for further layering.

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*The missing areas will be build-up with
incisal and tranparent porcelains.*

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* The firing parameters given above are guidelines, which must always be adjusted to suit the furnace used and the situation of the furnace. What is essential is getting the right firing result.

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*All tested materials conform to EN ISO 9693-1. The technical and
physical values quoted relate to samples produced in-house at our
own measuring instruments.*

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